# PRODUCT INFORMATION **Antibody Gel**

Catalog Number: AG-001-2

**Description:** Pure Goat Anti-Human IgG Gel, Specific Immobilized

Antibodies, 2mL

Lot Number:

**Bead Size:** 50-250 micron.

Spacer: None.

Linkage Between Spacer & Bead:

Imidoester.

Storage Buffer: 0.01M Phosphate - 0.15M NaCl, (PBS) pH 7.2-7.4. Contains

0.05% sodium azide as a preservative.

Storage: Store refrigerated at 5 - 8°C in Buffer. DO NOT FREEZE.

Stability: The material is stable for several years when stored

refrigerated with 0.05% sodium azide added as a preservative

Caution: Refer to the enclosed MSDS for information. The aluminum

seals have sharp edges and the vial itself may have cracks

which can cause lacerations. Use caution when opening the

vial.

# **EY** LABORATORIES, INC.

107 North Amphlett Blvd. San Mateo, CA 94401

650-342-3296 Tel: 650-342-2648 Orders: 1-800-821-0044 (Outside CA only)

## **GENERAL PROCEDURE Antibody Gel**

The following information is provided only for your convenience. The success of your experiments are not guaranteed by EY Laboratories, Inc.

# Use:

- **Procedure for** 1. The reaction may be carried out in a test tube (which will require centrifugation for washing steps) or in a small column (either a glass pipette or a plastic minicolumn.) Gels may be run at room temperature or in a cold room. Elevated temperatures should be avoided.
  - 2. Wash gel with 10 times the gel volume using 0.01M PBS.
  - 3. Apply sample and wash unbound material from column with Buffer. DO NOT OVERWASH!! Extensive washing may elute the protein to be purified if the binding constant is low.
  - 4. Elute bound material using 0.1M Glycine Buffer, pH 3.5. Collect small samples.
  - 5. Add 1.0M Tris-base in water to neutralize the acidic condition of the sample.

# Procedure for Gel Regeneration:

1. After elution, wash the gel with 10 times the gel volume using 1.0-1.4 M NaCl in distilled water. Re-equilibrate the gel by washing with 50 times the gel volume using Buffer (see reverse side). Store refrigerated with 0.05% sodium azide as a preservative. DO NOT store the column in the high salt concentration solution. DO NOT FREEZE.

**X** LABORATORIES, INC. 107 North Amphlett Blvd.

650-342-3296 Tel: 650-342-2648 Orders: 1-800-821-0044 San Mateo, CA 94401 (Outside CA only)

### MATERIAL SAFETY DATA SHEET

Effective Date: March 31, 2006 Revision 4 Page 1 of 2

#### PRODUCT IDENTIFICATION

Proteins, carbohydrates, and biotin immobilized on a support matrix of Name:

acrylamide or polygalactose.

ABP-01, A-1102 to A-9000, MB-1104 to MB-9000, PB-1104 to PB-Catalog Number (s):

9000, PG-001 to PG-7011, PB-01 to PB-05, CG-001 to CG-092, AG-

001 to AG-032, A-1001 to A-1004, CG-094 to CG-096.

Synonyms: Protein A, Avidin (egg white), D-Biotin, Lectins, Secondary

> Antibodies, Carbohydrates, Thyroglobulin, Fetuin, Hemoglobin, α-Lactalbumin, Porcine Stomach Mucin, Ovalbumin, Bovine Submaxillary Mucin, Transferrin, Myoglobin, Strept, Avidin, and 2-Iminobiotin immobilized on a polygalactose matrix or an acrylic (or

polyacrylaminde) matrix.

#### **EMERGENCY INFORMATION**

EY Laboratories, Inc. **EMERGENCY PHONE:** 107 North Amphlett Blvd. 650-342-3296 San Mateo, CA 94401

#### HAZARDOUS COMPONENTS

Specific protein(s) as listed on the vial label. The matrix itself is not known to be hazardous. The proteins are covalently attached to the beaded matrix and therefore present a hazard primarily through ingestion or injection. The biological activities of these chemicals will vary. It is possible that the immobilized material may leach off the beaded matrix during use. Care should be used when handling any of these reagents. All of these solutions contain at least 0.1%, but not greater than 1%, sodium azide as a preservative.

#### **HEALTH HAZARD INFORMATION**

EXPOSURE LIMITS: None established. The toxicological properties of these products have not

been thoroughly investigated. Care should be taken when handling any of

these materials.

**EFFECTS OF** No effects of overexposure have been documented. The individual proteins OVEREXPOSURE: and other ligands may cause allergic reactions in sensitive individuals. This

is a problem primarily with material that leaches from the column through

use. Local irritation is likely if eye contact occurs.

**ROUTES OF** Ingestion or injection of the beaded material are the primary routes of

EXPOSURE: exposure. Contact with the eyes may also present a hazard.

PHYSICAL CHARACTERISTICS

SOLUBILIZ

APPEARANCE: Solution containing a maximum of 50% (v/v) of beaded

matrix in buffer.

Not applicable.

**L** I LABORATORIES, INC.

Fax: 650-342-2648 107 North Amphlett Blvd. 1-800-821-0044 San Mateo, CA 94401 Orders:

(Outside CA only)

Tel: 650-342-3296

MSDS for Gels Continued - page 2 of 2.

FIRE AND EXPLOSION HAZARDS Not considered to be a fire hazard.

EXTINGUISHING MEDIA: Water spray or CO<sub>2</sub>.

SPECIAL FIRE FIGHTING PRECAUTIONS: None required.

NOTE: All solutions contain less than 1% sodium azide (w/v) as a preservative. Azide may react with lead and copper plumbing to form explosive metal azides. Flush with copious amounts of water when disposing material in the sink.

REACTIVITY DATA

STABILITY: Stable. Decomposition products are not known to be

hazardous

HAZARDOUS POLYMERIZATION: Will NOT occur.

INCOMPATIBILITY: None known. (Lead and copper may react with

SPILL / LEAK PROCEDURES

MATERIAL RELEASE / Avoid contact with liquid. Clean up spill with a paper towel soaked SPILL:

in household bleach. Do not allow solutions to dry on environmental surfaces. Wash affected area with detergent after the

area has been treated with bleach.

WASTE DISPOSAL: Incinerate, autoclave, or dispose of paper waste in accordance with

> all Local, State, and Federal regulations. Due to the small quantities of material involved these products are generally not considered to be environmental hazards. All of these proteins are

fully biodegradable.

#### **EMERGENCY FIRST AID PROCEDURES**

May be harmful if swallowed, injected, or allowed to contact the eyes. Wash contacted area with water for 15 minutes. If inhaled remove to fresh air, Report exposure to the appropriate safety official. Consult a physician if irritation occurs or if there is any indication of an allergic response, such as watering eyes, sneezing, or difficulty breathing. Any eye contact should be reported to a physician immediately.

#### SPECIAL HANDLING PRECAUTIONS

VENTILATION: No special ventilation is required but it is recommended to handle these

reagents in a fume hood when possible.

Required. Goggles or safety glasses with a side shield are Recommended EYE PROTECTION: RESPIRATORY Not required unless the formation of aerosols is likely. An approved PROTECTION:

respirator may be required for those individuals already known to be

sensitive to these materials.

**PROTECTIVE** Required when handling any of these materials.

GLOVES:

#### SPECIAL PRECAUTIONS

This material is for research and experimental application only. It is not intended for food, drug, household, agricultural, or cosmetic use. All materials should be handled only by technically qualified individuals experienced with working with potentially hazardous chemicals. The above information is correct to the best of our knowledge. The user should make independent decisions regarding completeness of the information, based on all sources available. EY Laboratories, Inc. shall not be held liable for any damage resulting from handling or contact with the above product.

**EY** LABORATORIES, INC.

107 North Amphlett Blvd. San Mateo, CA 94401

650-342-3296 Fax: 650-342-2648 1-800-821-0044 Orders:

(Outside CA only)