PRODUCT INFORMATION Colloidal Gold Conjugate

Catalog Number:	GAF-446-15 PLEASE FOLLO RECEIVED TO					
Description:	Colloidal Gold Conjugated Goat Purified Antibody To Mouse IgA, 1mL, 15nm Particle Size, 1mL	MATERIALS REQUIRED:	1.	InstantChel antigen (1µ		
Lot Number:			2.	80µl 1% B		
Expiration Date:			3.	20µl gold c		
Protein Concentration:				If less reage		
HAuCl₄ Concentration:	OD 520 = 4	INSTANTCHEK™	1.	Prepare a 1		
Buffer:	Stored in 0.002M Sodium Borate, pH 7.2 - 7.6	QC PROCEDURE FOR SPECIES IgG:	2.	Dilute 10µ dilute descr		
Storage:	DO NOT FREEZE ! Store at 5-8°C.		3.	Add the r device mem		
Other Data:	Absorbance peak at 520nm. CONTAINS 0.02% SODIUM AZIDE		4.	Wash the n		
	AS A PRESERVATIVE. Gently resuspend any sediment. If necessary, clarify by centrifugation at 400-500 x g for ten minutes before use.		5. ●	Result inter Dependin dots show		
Quality Control:	Procedure developed by EY Laboratories, Inc.			and passe		
References:	 Muller, C. et al. J. of Imm. Methods (1980) 37 :185-190. Roth, J. Techniques in Immunocytochemistry. Eds. Bullock and Petrusz Academic Press 1983, 217-284. 		•	No dots s have beer		
		NOTE:	Qu	estions or con		

INSTANTCHEK™ Q.C. PROCEDURE FOR GOLD CONJUGATES

PLEASE FOLLOW THE ACCOMPANYING PROCEDURE WHEN RECEIVED TO INSURE THAT THE MATERIAL IS ACTIVE.

MATERIALS REQUIRED:	1.	InstantChek Device impregnated with appropriate antigen $(1\mu g/spot)$ – included.				
	2.	80µl 1% BSA in PBS, pH 7.4.				
	3.	20µl gold conjugate.				
		If less reagent is used, please keep proportions constant.				
INSTANTCHEK™	1.	Prepare a 1% BSA in PBS solution.				
QC PROCEDURE FOR SPECIES IgG:	2.	Dilute 10µl Protein A Colloidal Gold with 90µl of dilute described in Step 1 and mix well.				
FOR 51 ECIE5 199.	3.	. Add the mix into the center of the InstantChek [™] device membrane and let absorb.				
	4.	Wash the membrane with $200\mu l$ of solution in Step 1.				
	5. ●	Result interpretation: Depending on the device, either one dot or all five dots show up. The colloidal gold conjugate is good and passes QC.				
	•	No dots show up. The reagent has problems and may have been destroyed during shipment.				
NOTE:	~	estions or comments? Please visit o://support.eylabs.net				

EY LABORATORIES, INC. 07 North Amphlett Blvd. San Mateo, CA 94401

Tel: 650-342-3296 Fax: 650-342-2648 Orders: 1-800-821-0044 (Outside CA only)



Tel: 650-342-3296 Fax: 650-342-2648 Orders: 1-800-821-0044 (Outside CA only)

MATERIAL SAFETY DATA SHEET

Effective Date: March 31, 2006 Revision 4 Page 1 of 2

PRODUCT IDENTIFICATION

 Name:
 Colloidal gold and colloidal silver labeled proteins, enzymes, and ligands.

 Catalog Number(s):
 G-2 to G-40, XGP-2, GP-01 to GP-8006, GAP-01, FGP-01, HGP-01, RGP-01, TGP-01, CCG-0001 to CCG-1018, GA-02, GAA-02, GAB-01 to GAB-02, FGA-02, HGA-02, GB-01 to GB-02, GE-01 to GE-03, GH-01 to GH-02, GM-01 to GM-2701, GAF-001 to GAF-2404, SA-02, SB-01, SH-01, SP-01 to SP-014.

 Formula:
 Complex polypeptides, enzymes, lectins, antibodies, and ligands coupled to colloidal gold or silver particles. Also, unconjugated colloidal gold particles.

Synonyms: Protein A, Horseradish Peroxidase, Strept. Avidin, D-Biotin, Purified Antibodies, Bovine Serum Albumin, Fetuin, Ovomucoid, RNase, DNase I, Alkaline Phosphatase, Protein G, Monoclonal Antibodies, Lectins, Neoglycoproteins, Adriamycin, and Neomycin coupled to colloidal gold particles or silver colloidal particles.

EMERGENCY INFORMATION

EY Laboratories, Inc. 107 North Amphlett Blvd. San Mateo, CA 94401

EMERGENCY PHONE: 650-342-3296

HAZARDOUS COMPONENTS

Specific protein or ligand as listed on the vial label. These solutions contain less than 0.1mg per ml. Biological activity of these proteins will vary. Although these materials are not generally considered to be hazardous they may cause allergic responses in sensitive individuals if inhaled or allowed to contact skin. Adriamycin and Neomycin are both used in cancer therapy and are cytoxic.

EXTREME CARE should be used when handling either of these two items. The colloidal gold and colloidal silver solutions are potentially caustic and will temporarily discolor the skin. Most solutions contain 0.02% sodium azide as a preservative.

HEALTH HAZARD INFORMATION

EXPOSURE LIMITS: None established. The toxicological properties of these products have not been thoroughly investigated. Care should be taken when handling any of these materials.
 EFFECTS OF OVEREXPOSURE: Any of these proteins may cause acute localized eye, skin, or mucous membrane irritation. Some sensitive individuals may develop a chronic allergic reaction with exposure.
 ROUTES OF EXPOSURE: Skin, eye, and mucous membrane contact. Care should be taken to avoid the formation of aerosols when handling any of these solutions.

PHYSICAL CHARACT RISTICS

APPEARANCE: The transmission of transmission of the transmission of transmission o

LABORATORIES, INC.

07 North Amphlett Blvd. San Mateo, CA 94401 Tel: 650-342-3296 Fax: 650-342-2648 Orders: 1-800-821-0044 (Outside CA only) MSDS for Colloidal Gold Labeled Proteins, Enzymes, and Ligands Continued - page 2 of 2.

FIRE AND EXPLOSION HAZARDS

EXTINGUISHING MEDIA: SPECIAL FIRE FIGHTING PRECAUTIONS: NOTE:

REACTIVITY DATA STABILITY:

HAZARDOUS POLYMERIZATION: INCOMPATIBILITY:

SPILL / LEAK PROCEDURES

MATERIAL RELEASE / SPILL:

WASTE DISPOSAL:

Water spray or CO₂. None required. Most solutions contain 0.05% sodium azide as a preservative. Azide may react with lead and copper plumbing to form explosive metal azides. Flush with copious amounts of water when disposing material in the sink.

Not considered to be a fire hazard.

Stable. Decomposition products are not known to be hazardous. Will NOT occur. None known. (Lead and copper may react with sodium azide).

Avoid contact with liquid. Clean up spill with a paper towel soaked in household bleach. Do not allow solutions to dry on environmental surfaces. Wash affected area with detergent after the area has been treated with bleach.

Incinerate, autoclave, or dispose of paper waste in accordance with all Local, State, and Federal regulations. Due to the small quantities of material involved these products are generally not considered to be environmental hazards. All of these proteins are fully biodegradable.

EMERGENCY FIRST AID PROCEDURES

May be harmful if swallowed, inhaled, or allowed to absorb through the skin. Wash contacted area with water for 15 minutes. If inhaled remove to fresh air. Report exposure to the appropriate safety official. The gold and silver sols may be caustic. Consult a physician if irritation occurs, if there is any indication of an allergic response such as watering eyes, sneezing, or difficulty breathing, or if eye contact occurs.

SPECIAL HANDLING PRECAUTIONS

VENTILATION:	No special ventilation is required but it is recommended to
EYE PROTECTION:	handle these reagents in a fume hood when possible. Safety goggles or safety glasses with side shields are recommended.
RESPIRATORY PROTECTION:	Recommended as a safety precaution. An approved respirator may be required for those individuals already known to be sensitive to these materials.
PROTECTIVE GLOVES:	Required when handling any of these materials.

SPECIAL PRECAUTIONS

This material is for research and experimental application only. It is not intended for food, drug, household, agricultural, or cosmetic use. All materials should be handled only by technically qualified individuals experienced with working with potentially hazardous chemicals. The above information is correct to the best of our knowledge. The user should make independent decisions regarding completeness of the information, based on all sources available. EY Laboratories, Inc. shall not be held liable for any damage resulting from handling or contact with the above product.



Tel: Fax:	650-342-3296 650-342-2648
Orders:	1-800-821-0044
	(Outside CA only)