

PRODUCT INFORMATION

Colloidal Gold Conjugate

Catalog Number: GAF-512-5

Description: Colloidal Gold Conjugated Rabbit Affinity Purified Antibody to Bovine IgG, Fc Specific, 5nm Particle Size, 1mL

Lot Number:

Expiration Date:

Protein Concentration: 10-20µg/ml

HAuCl₄ Concentration: OD 520 = 1.0

Buffer: Stored in 0.002M Sodium Borate, pH 7.2- 7.6

Storage: **DO NOT FREEZE** ! Store at 5-8° C.

Other Data: Absorbance peak at 520nm. CONTAINS 0.02% SODIUM AZIDE AS A PRESERVATIVE. Gently resuspend any sediment. If necessary, clarify by centrifugation at 400-500 x g for ten minutes before use.

Quality Control: Procedure developed by EY Laboratories, Inc.

References:

1. Muller, C. et al. J. of Imm. Methods (1980) **37** : 185-190.
2. Roth, J. Techniques in Immunocytochemistry. Eds. Bullock and Petrusz Academic Press 1983, 217-284.

INSTANTCHEK™ Q.C. PROCEDURE FOR GOLD CONJUGATES

PLEASE FOLLOW THE ACCOMPANYING PROCEDURE WHEN RECEIVED TO INSURE THAT THE MATERIAL IS ACTIVE.

MATERIALS REQUIRED:

1. InstantChek Device impregnated with appropriate antigen (1µg/spot) – included.
 2. 80µl 1% BSA in PBS, pH 7.4.
 3. 20µl gold conjugate.
- If less reagent is used, please keep proportions constant.

INSTANTCHEK™ QC PROCEDURE FOR SPECIES IgG:

1. Prepare a 1% BSA in PBS solution.
2. Dilute 10µl Protein A Colloidal Gold with 90µl of dilute described in Step 1 and mix well.
3. Add the mix into the center of the InstantChek™ device membrane and let absorb.
4. Wash the membrane with 200µl of solution in Step 1.
5. Result interpretation:
 - Depending on the device, either one dot or all five dots show up. The colloidal gold conjugate is good and passes QC.
 - No dots show up. The reagent has problems and may have been destroyed during shipment.

NOTE:

Questions or comments? Please visit <http://support.eylabs.net>

EY LABORATORIES, INC.
107 North Amphlett Blvd.
San Mateo, CA 94401

Tel: 650-342-3296
Fax: 650-342-2648
Orders: 1-800-821-0044
(Outside CA only)

EY LABORATORIES, INC.
107 North Amphlett Blvd.
San Mateo, CA 94401

Tel: 650-342-3296
Fax: 650-342-2648
Orders: 1-800-821-0044
(Outside CA only)

Sample Only

MATERIAL SAFETY DATA SHEET

Effective Date: March 31, 2006

Revision 4

Page 1 of 2

PRODUCT IDENTIFICATION

Name: Colloidal gold and colloidal silver labeled proteins, enzymes, and ligands.

Catalog Number(s): G-2 to G-40, XGP-2, GP-01 to GP-8006, GAP-01, FGP-01, HGP-01, RGP-01, TGP-01, CCG-0001 to CCG-1018, GA-02, GAA-02, GAB-01 to GAB-02, FGA-02, HGA-02, GB-01 to GB-02, GE-01 to GE-03, GH-01 to GH-02, GM-01 to GM-2701, GAF-001 to GAF-2404, SA-02, SB-01, SH-01, SP-01 to SP-014.

Formula: Complex polypeptides, enzymes, lectins, antibodies, and ligands coupled to colloidal gold or silver particles. Also, unconjugated colloidal gold particles.

Synonyms: Protein A, Horseradish Peroxidase, Strept. Avidin, D-Biotin, Purified Antibodies, Bovine Serum Albumin, Fetuin, Ovomucoid, RNase, DNase I, Alkaline Phosphatase, Protein G, Monoclonal Antibodies, Lectins, Neoglycoproteins, Adriamycin, and Neomycin coupled to colloidal gold particles or silver colloidal particles.

EMERGENCY INFORMATION

EY Laboratories, Inc.
107 North Amphlett Blvd.
San Mateo, CA 94401

EMERGENCY PHONE:**650-342-3296****HAZARDOUS COMPONENTS**

Specific protein or ligand as listed on the vial label. These solutions contain less than 0.1mg per ml. Biological activity of these proteins will vary. Although these materials are not generally considered to be hazardous they may cause allergic responses in sensitive individuals if inhaled or allowed to contact skin. Adriamycin and Neomycin are both used in cancer therapy and are cytotoxic.

EXTREME CARE should be used when handling either of these two items. The colloidal gold and colloidal silver solutions are potentially caustic and will temporarily discolor the skin. Most solutions contain 0.02% sodium azide as a preservative.

HEALTH HAZARD INFORMATION

EXPOSURE LIMITS: None established. The toxicological properties of these products have not been thoroughly investigated. Care should be taken when handling any of these materials.

EFFECTS OF OVEREXPOSURE: Any of these proteins may cause acute localized eye, skin, or mucous membrane irritation. Some sensitive individuals may develop a chronic allergic reaction with exposure.

ROUTES OF EXPOSURE: Skin, eye, and mucous membrane contact. Care should be taken to avoid the formation of aerosols when handling any of these solutions.

PHYSICAL CHARACTERISTICS

APPEARANCE: Light burgundy to purple liquid. 2nm - pale yellowish-brown liquid.

SOLUBILITY: All liquids are completely miscible in water and biological buffers.

FIRE AND EXPLOSION HAZARDS

Not considered to be a fire hazard.

EXTINGUISHING MEDIA:Water spray or CO₂.**SPECIAL FIRE FIGHTING PRECAUTIONS:**

None required.

NOTE:

Most solutions contain 0.05% sodium azide as a preservative. Azide may react with lead and copper plumbing to form explosive metal azides. Flush with copious amounts of water when disposing material in the sink.

REACTIVITY DATA**STABILITY:**

Stable. Decomposition products are not known to be hazardous.

HAZARDOUS POLYMERIZATION:

Will NOT occur.

INCOMPATIBILITY:

None known. (Lead and copper may react with sodium azide).

SPILL / LEAK PROCEDURES**MATERIAL RELEASE / SPILL:**

Avoid contact with liquid. Clean up spill with a paper towel soaked in household bleach. Do not allow solutions to dry on environmental surfaces. Wash affected area with detergent after the area has been treated with bleach.

WASTE DISPOSAL:

Incinerate, autoclave, or dispose of paper waste in accordance with all Local, State, and Federal regulations. Due to the small quantities of material involved these products are generally not considered to be environmental hazards. All of these proteins are fully biodegradable.

EMERGENCY FIRST AID PROCEDURES

May be harmful if swallowed, inhaled, or allowed to absorb through the skin. Wash contacted area with water for 15 minutes. If inhaled remove to fresh air. Report exposure to the appropriate safety official. The gold and silver sols may be caustic. Consult a physician if irritation occurs, if there is any indication of an allergic response such as watering eyes, sneezing, or difficulty breathing, or if eye contact occurs.

SPECIAL HANDLING PRECAUTIONS**VENTILATION:**

No special ventilation is required but it is recommended to handle these reagents in a fume hood when possible.

EYE PROTECTION:

Safety goggles or safety glasses with side shields are recommended.

RESPIRATORY PROTECTION:

Recommended as a safety precaution. An approved respirator may be required for those individuals already known to be sensitive to these materials.

PROTECTIVE GLOVES:

Required when handling any of these materials.

SPECIAL PRECAUTIONS

This material is for research and experimental application only. It is not intended for food, drug, household, agricultural, or cosmetic use. All materials should be handled only by technically qualified individuals experienced with working with potentially hazardous chemicals. The above information is correct to the best of our knowledge. The user should make independent decisions regarding completeness of the information, based on all sources available. EY Laboratories, Inc. shall not be held liable for any damage resulting from handling or contact with the above product.

EY LABORATORIES, INC.

107 North Amphlett Blvd.
San Mateo, CA 94401

Tel: 650-342-3296
Fax: 650-342-2648
Orders: 1-800-821-0044
(Outside CA only)

EY LABORATORIES, INC.

107 North Amphlett Blvd.
San Mateo, CA 94401

Tel: 650-342-3296
Fax: 650-342-2648
Orders: 1-800-821-0044
(Outside CA only)

Sample Only