PRODUCT INFORMATION Alkaline Phosphatase Labeled Lectins

Catalog Number: LA-5601-1

Description: Pure *Trichosanthes kirilowii* lectin (TKA) from China gourd,

tianhuafen, Alkaline Phosphatase conjugated.

Lot Number:

Protein Concentration:

1 mg purified TKA Alkaline Phosphatase / 1 ml Buffer.

Carbohydrate Specificity:

Galactose

Inhibitory Carbohydrate:

Lactose $> \beta$ -Galactose $> \alpha$ -Galactose

Activity: Less than 5 μ g/ml will agglutinate human type O erythrocytes. Less

than 0.5 µg/ml will agglutinate neuraminidase treated erythrocytes.

Buffer: 0.01M Phosphate - 0.15M NaCl, pH 7.2 - 7.4.

Chemical Used for Conjugation:

Alkaline Phosphatase.

Storage: Store liquid refrigerated at 5-8°C in aliquots. DO NOT FREEZE!

(20-50% Glycerol has been added to prevent freezing)

Stability: The liquid material is stable for at least 1 year when stored refrigerated

in aliquots with 0.05% sodium azide added as a preservative.

Caution: Refer to the enclosed MSDS for information regarding lectins. The

aluminum seals have sharp edges and the vial itself may have cracks which can cause lacerations. Use caution when opening the vial.

Procedure for Use: See reverse side.

References:

Yeung, et al. (1980) Proc. Asian Symp. Medicinal Plants and Spices.

Tel:

650-342-3296

4th, September, Bangkok, Thailand.

Schalla, W. O., et al (1985) J. Clin. Micro. 22: 379-382.

LABORATORIES, INC.

107 North Amphlett Blvd.
San Mateo, CA 94401

Fax: 650-342-2648
Orders: 1-800-821-0044
(Outside CA only)

PRODUCT INFORMATION Alkaline Phosphatase Enzyme Activity Assay

Chemical Principle:

Orthophosphoric Monoester + H₂0 → Alcohol + H₃PO₄

Assay Reagents: BUFFER: 0.1 M Tris buffer, pH 8.2.

ENZYME: Dilute with 0.1 M Tris Buffer.

Acceptable dilution: 5-20 µg/ml.

En

SUBSTRATE: 0.001 M p-nitrophenylphosphate (P-NPP).

Procedure: 1. Add 2.9 ml substrate to Reaction test tube and 2.9 ml to

Control test tube.

2. At time = 0, add 100 μl of diluted ENZYME to Reaction tube

and 100µl Tris to Control tube. Mix thoroughly.

3. Measure and record optical density at 410 nm OD(410) every 15 seconds for 3 minutes, or take end point reading after 3 minutes by stopping reaction with 100µl of 5.0 M NaOH.

4. Use the OD(410) measurement to determine the rate of change

in absorbance per minute.

Enzyme Activity Calculations:

One unit of activity is the amount of enzyme to decompose 1 µmole of P-NPP/minute at 25°C. 1.62 X 10⁴ cm⁻¹ is the molar absorbance of P-NPP.

accordance of F 1411.

OD(410) / min =

OD(410) / 3min - OD(410) Control / 3 minutes

3 minutes

mg enzyme / ml reaction mixture = [enzyme dilution]

,

units / mg = $\frac{OD(410) / min}{1.62 \times 10^4 ml reaction mixture}$

LABORATORIES, INC. 107 North Amphlett Blvd. San Mateo, CA 94401

Tel: 650-342-3296 Fax: 650-342-2648 Orders: 1-800-821-0044 (Outside CA only)

MATERIAL SAFETY DATA SHEET

Effective Date: March 31, 2006 Revision 4 Page 1 of 2

PRODUCT IDENTIFICATION

Purified proteins or biotin labeled with Horseradish Peroxidase or Alkaline Name:

Phosphatase.

Catalog HP-02, BA-104, BA-105, BA-108, BA-109, H-1102 to H-9000, LA-1104 to LA-9000, PA-2100 to PA-2701, AA-2100 to AA-2701, HAF-001 to HAF-2354, AAF-Number (s):

001 to AAF-2354, HA-01 to HA-013, AA-01 to AA-013, HAL-1104 to HAL-4701.

AAL-1104 to AAL-4701.

Protein A, Avidin (egg white), Biotin, Lectins, Secondary Antibodies labeled with Synonyms:

Horseradish Peroxidase or Alkaline Phosphatase.

EMERGENCY INFORMATION

EY Laboratories, Inc. **EMERGENCY PHONE:** 107 North Amphlett Blvd. 650 342 3296 San Mateo, CÂ 94401

HAZARDOUS COMPONENTS

Specific protein(s) as listed on the vial label. Solutions are at a concentration generally greater than 0.5mg protein/ ml. Biological activity of these labeled proteins will vary. Horseradish Peroxidase and Alkaline Phosphatase are both potent enzymes which may be harmful if ingested, inhaled, or allowed to absorb through the skin. Both enzymes are known to cause allergic responses in sensitive individuals

HEALTH HAZARD INFORMATION

EXPOSURE LIMITS. None established. The toxicological properties of these products have not

been thoroughly investigated. Care should be taken when handling any

of these materials.

EFFECTS OF May causes localized eye, skin, or mucous membrane irritation. Some OVEREXPOSURE:

sensitive individuals may develop a chronic allergic reaction with

ROUTES OF Inhalation of powders and skin contact with liquids are the primary **EXPOSURE:**

routes of exposure. Care should be taken to avoid the formation of

aerosols when handling any of the solutions.

PHYSICAL CHARACTERISTICS

APPEARANCE: Powders are a light brown. Solutions will be light to dark brown. SOLUBILITY: Powders are completely soluble in many biological buffers and water. All liquids are completely miscible in water and biological buffers.

FIRE AND EXPLOSION HAZARDS

Not considered to be a fire hazard. EXTINGUISHING MEDIA: Water spray or CO2. SPECIAL FIRE FIGHTING None required. PRECAUTIONS.

NOTE:

Alkaline Phosphatase conjugates contain less than 0.05% sodium azide as a preservative. Azide may react with lead and copper plumbing to form explosive metal azides. Flush with copious amounts of water when disposing material in the sink.

LABORATORIES, INC.

Tel: 650-342-3296 Fax: 650-342-2648 North Amphlett Blvd. 1-800-821-0044 San Mateo, CA 94401 Orders: (Outside CA only) MSDS for Horseradish or Alkaline Phosphatase Labeled Proteins & Biotin Continued - page 2 of 2.

REACTIVITY DATA

STABILITY: Stable. The nature of any decomposition products are

not known. They are not believed to be hazardous.

HAZARDOUS POLYMERIZATION: Will NOT occur.

INCOMPATIBILITY. None known. (Lead and copper may react with

sodium azide).

SPILL / LEAK PROCEDURES

MATERIAL RELEASE /

SPILL:

Avoid contact with powder or liquid. Clean up spill with a paper towel soaked in household bleach. Do not allow solutions to dry on environmental surfaces. Wash affected area with detergent after the area

has been treated with bleach.

Incinerate, autoclave, or dispose of paper waste in accordance with all WASTE DISPOSAL:

Local, State, and Federal regulations. Due to the small quantities of material involved these products are generally not considered to be environmental hazards. All of these proteins are fully biodegradable.

EMERGENCY FIRST AID PROCEDURES

May be harmful if swallowed, inhaled, or allowed to absorb through the skin. Wash contacted area with water for 15 minutes. If inhaled remove to fresh air. Report exposure to the appropriate safety official. Consult a physician if irritation occurs or if there is any indication of an allergic response, such as watering eyes, sneezing, or difficulty breathing. Any eye contact should be reported to a physician immediately

SPECIAL HANDLING PRECAUTIONS

VENTILATION: No special ventilation is required but it is recommended to handle these reagents in a fume hood when possible.

EYE PROTECTION: Required, Goggles or safety glasses with a side shield are

recommended.

RESPIRATORY PROTECTION: Recommended as a safety precaution, specifically when

> working with powders. An approved respirator may be required for those individuals already known to be

sensitive to these materials.

PROTECTIVE GLOVES: Required when handling any of these materials.

SPECIAL PRECAUTIONS

This material is for research and experimental application only. It is not intended for food, drug, household, agricultural, or cosmetic use. All materials should be handled only by technically qualified individuals experienced with working with potentially hazardous chemicals. The above information is correct to the best of our knowledge. The user should make independent decisions regarding completeness of the information, based on all sources available. EY Laboratories, Inc. shall not be held liable for any damage resulting from handling or contact with the above product.

LABORATORIES, INC.

107 North Amphlett Blvd. San Mateo, CA 94401

Tel: 650-342-3296 Fax: 650-342-2648 1-800-821-0044 Orders: (Outside CA only)